# INSTRUCTIONS FOR USE

---

## OSSURE™ Local Osteo-Enhancement Procedure (LOEP) Kit

**REF** AG11040302C

**Basic UDI-DI:** 0850029099AG11040302CFF

The following languages are included in this packet:

English (EN)	Deutsch (DE)	Français (FR)	Italiano (IT)
Español (ES)	Dansk (DK)	Nederlands (NL)	

This IFU is available online at: <https://www.agnovos.com/IFU>

For additional information and translations please contact the manufacturer or local distributor.

Customer Service Germany  +49 618 282 248 00  [customerservice@agnovos.com](mailto:customerservice@agnovos.com)

Customer Service Switzerland  +41 43 547 01 10  [customerservice@agnovos.com](mailto:customerservice@agnovos.com)

**EC** **REP** Emergo Europe  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

AgNovos Healthcare Germany GmbH  
Eschersheimer Landstrasse  
14 60322 Frankfurt am Main  
Germany

**CH** **REP** MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

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**CE**  
2797

**RxONLY**

The CE-Marking of Conformity is applied per catalog number and appears on the outer label, if applicable.

*Attention Operating Surgeon*

## **IMPORTANT MEDICAL INFORMATION**

**AGNOVOS HEALTHCARE  
OSSURE LOEP Kit**

---

**REF AG11040302C**

---

### **CONTENT**

DEFINITIONS


GENERAL PRODUCT INFORMATION

- A. INTENDED PURPOSE
- B. SUMMARY OF CLINICAL BENEFITS
- C. INDICATIONS
- D. CONTRAINDICATIONS
- E. INTENDED PATIENT POPULATION
- F. INTENDED USERS
- G. PERFORMANCE CHARACTERISTICS
- H. POTENTIAL COMPLICATIONS
- I. POTENTIAL EXPOSURE TO MATERIALS
- J. WARNING AND PRECAUTIONS
- K. ADVERSE EFFECTS
- L. HANDLING AND STERILIZATION
- M. STORAGE CONDITIONS
- N. OTHER GENERAL PURPOSE EQUIPMENT
- O. DIRECTIONS FOR USE / MIXING INSTRUCTIONS
- P. PATIENT INFORMATION LEAFLET & IMPLANT CARD
- Q. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)
- R. REPORTING SERIOUS INCIDENTS

# DEFINITIONS

Symbols and abbreviations may be used on the package label. The following table provides the definition of these symbols and abbreviations.

**Table 1: Definitions of Symbols and Abbreviations**

SYMBOL	DEFINITION
	CE Mark
	Medical device
	Batch code
	Catalog number
	Unique Device Identifier
	Do not reuse
	Do not re-sterilize
	Caution
	Consult instructions for use or consult electronic instructions for use
	Use by
	Storage temperature limitation
	Keep dry
	Keep away from sunlight
	Date of manufacture
	Name and address of manufacturer
	Sterilized using irradiation
	For prescription use only
	Contains hazardous substances
	Do not use if packaging is damaged and consult instructions for use
	Authorized Representative in the European Community / European Union
	Authorized Representative in Switzerland
	Importer

SYMBOL	DEFINITION
	Single sterile barrier system with protective packaging inside
	Single sterile barrier system with protective packaging outside
	Patient identification
	Patient information website
	Health care center or doctor
	Date
	MRI Safe
	Recycling symbol for High-Density Polyethylene (HDPE)
	Peel here
	Continued on reverse side
	Less than 5 minutes
	Question

## **GENERAL PRODUCT INFORMATION**

The OSSURE LOEP Kit is a sterile, single-use medical device kit containing the OSSURE implant material, the instruments needed to mix the implant material, and instruments helpful to access and prepare the implantation site and deliver the implant material. The OSSURE implant material is a calcium salt-based bone void filler consisting of a powder component and an aqueous mixing solution that are combined to form an injectable paste that hardens and cures in situ. Once implanted, the various calcium sulfate and calcium phosphate phases provide a resorbable scaffold for bone formation. Implant material resorption occurs through chemical dissolution and cell-mediated resorption. The resorption process is coupled with the formation of new bone. The resorption and bone formation process is expected to occur over a period of 6 to 24 months.

The components of the OSSURE LOEP Kit are packaged in two high density polyethylene sealed thermoformed trays. The upper tray, referred to as the Instrument Tray, includes the instruments that may be used to access and prepare the implantation site. The lower tray, referred to as the "Mixing Tray", includes the OSSURE implant material components, mixing instruments, and instruments that may be used to deliver the implant material. The "Instrument Tray" and the "Mixing Tray" are provided within a third, larger high- density polyethylene sealed thermoformed "Outer Tray", the primary sterile barrier, terminally sterilized using gamma radiation.

### **A. INTENDED PURPOSE**

The OSSURE resorbable implant material contained in the OSSURE LOEP Kit is intended to form new bone in voids not intrinsic to the structural stability of the skeletal system. The OSSURE LOEP Kit includes the instrumentation to mix the implant material, to access and prepare the implantation site, and to deliver the implant material.

### **B. SUMMARY OF CLINICAL BENEFITS**

- Forms new bone in skeletal voids - replacement of bone lost due to surgical procedures, trauma, or disease.
- Injectable material allows for use in minimally invasive procedures - procedure is less invasive, minimizing patient time to recovery.
- Hardens in situ – material provides biomechanical support while new bone is being formed.
- Material is synthetic - eliminates the need for autologous graft harvest and potential complications (e.g., pain, discomfort, bruising).
- Radiopaque - allows for fluoroscopic monitoring during injection.

### **C. INDICATIONS**

The OSSURE resorbable implant material contained in the OSSURE LOEP Kit is intended to form new bone in voids in the skeletal system. It is intended to be injected into surgically prepared sites in the pelvis and extremities, where it cures in situ, and is resorbed and replaced with new bone. OSSURE is appropriate for use in voids not intrinsic to the structural stability of the skeletal system.

## **D. CONTRAINDICATIONS**

The OSSURE resorbable implant material is contraindicated where the device is used in articulating surfaces, or for structural support in lieu of hardware in load-bearing bone.

Conditions representing relative contraindications include:

- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative bone disease that compromises the body's ability to replace the implant material with bone
- Closed voids/gaps that are structurally capable of over-pressurization during injection
- Pregnancy
- Uncooperative patients who will not or cannot follow postoperative instructions, including individuals who abuse drugs and/or alcohol
- Pre-existing calcium metabolism disorder (e.g., hypercalcaemia)
- Renal compromised patients
- Patients with a history of or active tuberculous spondylitis

The OSSURE LOEP Kit is contraindicated for use in spinal, cranial, and craniomaxillofacial reconstruction applications.

## **E. INTENDED PATIENT POPULATION**

The intended target population for the OSSURE LOEP Kit is adults with bone voids resulting from surgical procedures, trauma, or disease.

## **F. INTENDED USERS**

The OSSURE LOEP Kit is intended for use by medical practitioners qualified to perform orthopedic procedures.

## **G. PERFORMANCE CHARACTERISTICS**

Once mixed, the OSSURE implant material has a recommended working time of 5-7 minutes, not to exceed 8 minutes. The implant material is radiopaque. Upon implantation, the material hardens in situ in approximately 20 minutes, reaches more than 85% of its ultimate strength within 1 hour and is replaced with new bone over time. The expected resorption time is 6-24 months.

## **H. POTENTIAL COMPLICATIONS**

Proper surgical procedures and techniques are the responsibility of the medical professional. Each surgeon must evaluate the appropriateness of the procedure used based on personal medical training and experience.

## I. Potential Exposure to Materials

### OSSURE Implant Material

Each OSSURE LOEP Kit contains the powder and liquid components needed to produce up to 25cc of OSSURE Implant Material. The powder components consist of calcium salts that are packaged in a sealed, polycarbonate dual powder container. The neutralized glycolic acid liquid component is packaged in a capped glass vial. The components are as follows:

Material	Substance
Powder Component (56g)	Calcium Sulfate Salts (74.9%)
	Beta-Tricalcium Phosphate ( $\beta$ -TCP) (18.3%)
	Monocalcium Phosphate Monohydrate (MCPM) (6.7%)
	Sucrose (0.1%)
Liquid Component (13g)	USP Water (93.5%)
	Glycolic Acid (4.3%)
	Sodium Hydroxide (2.2%)

Once the powder and liquid components are mixed, the calcium sulfate salts hydrate with water to form calcium sulfate dihydrate. At the same time, the MCPM and  $\beta$ -TCP dissolve in the liquid and react to form brushite.

During these reactions, a small amount of phosphoric acid is produced, which is quickly consumed. Any remaining water and glycolic acid are incorporated into the crystal structure of the set implant material.

After mixing and setting, the OSSURE Implant Material is a tri-phasic, calcium-based cement consisting of the following:

- 75% Calcium Sulfate (52g)
- 9% Beta Tricalcium Phosphate ( $\beta$ -TCP) (6g)
- 16% Brushite (11g)

Upon implantation, the material sets in approximately 20 minutes and reaches more than 85% of its ultimate strength within 1 hour.

### Stainless Steel Surgical Instruments

The surgical instruments provided in the OSSURE LOEP kit (i.e. the Suction Irrigator, Injection Cannula and Serrated Tissue Protector) are manufactured from 304 stainless steel and contain trace amounts [ $>0.1\%$  (wt)] of nickel and cobalt as an alloying element.

## J. WARNINGS AND PRECAUTIONS

As with any surgical procedure, care should be exercised in treating individuals with preexisting conditions that may affect the success of the surgical procedure. This includes individuals with bleeding disorders of any etiology, cardiovascular disease, long-term steroidal therapy, immunosuppressive therapy, or high dosage radiation therapy.

The procedure to implant the OSSURE implant material may be conducted under fluoroscopic guidance to monitor implantation resulting in radiation exposure to the patient.

Avoid overfilling or pressurizing the implantation site to prevent embolization of fat and/or implant material.

Use the OSSURE implant material as supplied and according to the Handling and Use information provided. All powder (supplied) and all solution (supplied) must be used when mixing the implant material.

**Warning: Do not use kit if the sterile packaging is damaged or unintentionally opened before use.**

Do not add any additional substances to the OSSURE implant material. Performance characteristics of this product are based on the formulation provided. Using alternate mixing solutions and/or adding other substances to the mixture may alter the curing and performance of this product.

The OSSURE implant material is non-metallic. It is not expected to cause interference with typical medical examinations or be impacted by environmental conditions or other external sources.

The OSSURE implant material is MRI safe.

### **Intra-Operative Precautions**

- Inspect instruments prior to use for damage during shipment or storage or any out-of-box defects. Do not use any instruments that appear damaged and replace with an instrument from a new OSSURE LOEP Kit. If possible, retain and photograph the damaged instrument to assist with AgNovos Healthcare's investigation of the deficiency.
- Inspect instruments immediately upon removal from the patient for any signs of damage. If an instrument appears damaged, discontinue use. If necessary, replace with an instrument from a new OSSURE LOEP Kit. If possible, retain and photograph the damaged instrument to assist with AgNovos Healthcare's investigation of the deficiency.
- The instruments provided in the kit are intended for single use. Should any procedural steps need to be repeated, used instruments should be replaced with instruments from new OSSURE LOEP Kit.
- If material extravasates laterally into the soft tissue, use a moistened sponge and pickups to remove implant material.

## **K. ADVERSE EFFECTS**

Possible adverse effects include but are not limited to:

- Wound complications including bruising, hematoma, site drainage, bone fracture, infection, and other complications that are possible with any surgery
- Fracture or extrusion of the OSSURE implant material, with or without particulate debris generation
- Deformity of the bone at the site
- Incomplete, or lack of, osseous ingrowth into the implantation site
- Transient hypercalcaemia
- Potential to over-pressurize the implantation site which could result in embolization of fat and/or implant material

The occurrence of any of these effects may require medical and/or surgical intervention.

## **Residual Risks & Undesirable Side Effects**

Examples of residual risks, specific to the implantation of bone void fillers and bone cements, that should be considered include:

- Osteolysis
- Wound infection
- Fat emboli and sequelae
- Venous thrombosis and sequelae
- Material extravasation into vessels / emboli and sequelae
- Tissue irritation
- Treated bone fracture
- Injury to blood vessels
- Pain
- Wound dehiscence

## **L. HANDLING AND STERILIZATION**

The OSSURE implant material is provided sterile via gamma radiation and should be considered sterile unless the primary sterile barrier packaging has been opened or damaged. The kits are for single patient use and should never be re-sterilized or reused. Reuse of these devices may potentially result in serious patient harm. Examples of hazards related to the reuse of these devices include but are not limited to significant degradation in device performance, cross infection, and contamination.

If the sterile packaging has been compromised, the kit materials shall be considered contaminated and shall be disposed of in accordance with local requirements. No additional special handling or disposal requirements apply.

**CAUTION: Surgical instrumentation contained in the kit may contain sharp edges. Care should be taken when handling. Dispose of properly.**

## **M. STORAGE CONDITIONS**

All kits must be stored in a clean, dry environment and be protected from sunlight and extremes in temperature. The OSSURE LOEP Kit should be stored at 15-30°C or 59-86°F.

## **N. OTHER GENERAL PURPOSE EQUIPMENT**

Drilling instruments for orthopedic procedures are required to access the implant site and must be used strictly according to their intended purpose. The following instruments are also required:

- 2.5mm-compatible pin collet chuck (not included)
- 6.35mm chuck (not included)


## **O. DIRECTIONS FOR USE / MIXING INSTRUCTIONS**

The OSSURE LOEP Kit contains the components and instruments used to mix and inject the resultant implant material.

### **MIX OSSURE IMPLANT MATERIAL**

- DO NOT mix OSSURE implant material until the surgeon has debrided and irrigated the enhancement site.
- Follow the mixing instructions exactly.
- The curing process begins when the powder components are combined; do not pre-mix the powders.
- Mixing and syringe filling should be completed in less than 90 seconds.
- Injection of the material should be completed within 5-7 minutes from when the two powder components are mixed.

1. Grasp the mixer top (**Figures A, B**), and turn it counter-clockwise to loosen; remove it fully from mixer. Carefully set mixer body and mixer top aside being careful not to soil either part.

2. Stand powder container (**Figure C**) upright on flat surface; remove lid. Place mixer on top of powder container. Twist mixer clockwise to engage threads and dock the components.

**FIGURE C**

**NOTE: Do not start mixing until surgeon has completed debridement/suction and irrigation of the enhancement site and the second timeout has been completed.**

**NOTE: The team member mixing should announce commencement of mixing.**

3. Turn over container and mixer to add all powder from container into mixer body; (**Figure D**) execute one more full inversion cycle. When complete, the mixer should be below the powder container.

**NOTE: Tap the container as needed and ensure that all powder has transferred to the mixer.**

Disconnect powder container from mixer and discard.

**FIGURE D**

**NOTE: Setting of OSSURE begins when the two powder components have been mixed. Ideally, injection should be completed within 5-7 minutes.**

4. Carefully pour all mixing solution into the mixer body (**Figure E**).

NOTE: DO NOT add any additional substances to OSSURE implant material. Performance characteristics of this product are based on the formulation provided. Using alternate mixing solutions and/or adding other substances to the mixture may alter the curing and performance of this product.

5. Reassemble the mixer and turn handle vigorously in a single direction for **30 seconds (Figure F)**; during that time execute **at least 50 full revolutions**.

6. Remove and retain red tab (**Figure G**).

7. Remove mixer top from mixer body; if implant material is stuck to mixer top blade, use red tab from mixer body and scrape as much material as possible into mixer body (**Figure H**).

**FIGURE H**

8. Turn mixer body so that it is parallel to the floor and rotate in hand until syringe port is facing you. Grasp the base of the mixer and twist it towards you until the mixer floor is fully lowered (base of the mixer body no longer twists and the syringe port is fully exposed) (**Figure I**).

**FIGURE I**


9. Locate the threaded extruder and insert it into the mixer body (**Figure J**), pushing it until completely seated into mixer.

**FIGURE J**

10. Simultaneously push and twist the threaded extruder clockwise (**Figure K**) until there is no air space between it and OSSURE implant material; continue twisting until the material reaches just beyond the end of the syringe port (**Figure L**).

**NOTE: Be careful not to advance the threaded extruder so far that the implant material escapes fully from the syringe port resulting in loss of material.**

**FIGURE K**


**FIGURE L**

11. Remove and retain white cap from syringe (**Figure M**). Insert the end of the syringe barrel onto the syringe port, accomplishing a secure friction fit (**Figure N**).

**FIGURE M**

**FIGURE N**

12. Advance the threaded extruder until approximately 12 mL of implant material is in the syringe (**Figure O**).

**FIGURE O**

**NOTE: Ensure syringe remains connected to mixer while advancing the extruder without overfilling the syringe, potentially resulting in the loss of material.**

13. Re-attach the white cap to syringe (**Figure P**) and then attach the injection cannula (**Figure Q**) to the syringe and depress the plunger slowly to fill the cannula with OSSURE being careful to avoid loss of material; hand syringe to the surgeon when prompted.

**FIGURE P**

**FIGURE Q**

**NOTE: When attaching the cannula to the syringe, avoid over-tightening as this could lead to breaking the injection cannula.**

14. Immediately repeat steps 11-13 filling the second syringe with the remaining implant material while the surgeon uses the first syringe; total volume added to syringes should be between 20-25 mL.

## **P. PATIENT INFORMATION LEAFLET & IMPLANT CARD**

Each patient treated with OSSURE LOEP should be provided with the Patient Information Leaflet and Implant Card that are included in the kit. The treating physician must enter the following information on the card before providing to the patient:

- Patient Name or ID
- Date of Implantation
- Name and Address of the Healthcare Institution or Provider

## **Q. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)**

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

- URL to the Eudamed public website: <https://ec.europa.eu/tools/eudamed>.
- Basic UDI-DI: 0850029099AG11040302CFF

You can also request a copy of the SSCP by contacting AgNovos Customer Service at the number or email address below.

## **R. REPORTING INCIDENTS**

Serious incidents related to use of the OSSURE LOEP Kit should be reported to AgNovos Customer Service and to the Competent Authority of the country where the device is used.

**AgNovos Customer Service Europe**  
**+49 618 282 248 00**  
**customerservice@agnovos.com**

## OSSURE™ Local Osteo-Enhancement Procedure (LOEP) Set

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